



Titrette®
Bottle-top burette, class A precision

Bottle-top burette Titrette®, DE-M, 50ml, Standard
Minimum order quantity / quantity interval: 1 piece

Cat. No.: 4760161
705,60 EUR

Recommended sales price without VAT (PL 2018)
per 1 on orders of 1.

Product description

Titrette® - the first bottle-top burette to satisfy Class A error limits for glass burettes!

It is also remarkable for its easy dropwise titration, compact design, ultra-high precision, replaceable piston/cylinder assembly and an optional PC interface. With the Titrette® bottle-top burette, you can titrate quickly and reliably with the highest precision, even in close quarters, with no power hookup needed - in the lab, in production, or in the field.

Additional features:

- The control elements of the Titrette® bottle-top burette have an intuitive layout.
- Separate buttons for On/Off and Pause; CLEAR button allows user to reset the display and select functions.
- Easy-grip hand wheel action and smooth precision gears for fast or drop-wise titration make handling more reliable and extremely simple.

Items supplied

Titrette® bottle-top burette, DE-M marking, supplied with performance certificate, telescoping filling tube (170 - 330 mm), recirculation tube, 2 batteries (AAA/UM4/LR03), 3 PP bottle adapters (GL 45/32, GL 45/S 40, GL 32/NS 29/32), 2 amber colored light shield inspection windows. Additionally included with RS 232 interface: 2 m interface cable (Sub-D plug connector, 9-pin), one CD (driver software and open RS232 communication protocol). The CD also includes an example application in XLS-file format, as well as a special operating manual.

Specification

Capacity	ml	50
A (%)	≤ ± %	0,06
A (µl)	≤ ± µl	30
CV (%)	≤ %	0,02
CV (µl)	≤ µl	10
version		Standard